

M-Commander InVision

Seamless Transparency




M-Commander InVision Analysis Software
for hygiene, logistics and technology



Keep track to keep safe

Discover the benefits of having a seamless, centralised and permanent overview of all the relevant data at all levels. Lightning-fast data transfers across large distances ensure the efficient management of hygiene, logistics and service – putting things firmly under your control.

M-Commander
InVision



Hygiene management and cycle management

Efficient information management

All the cycle data relating to your MEIKO cleaning and disinfection machines is continuously documented. You can compile this information into a handy log report whenever you like, and then print it or store it on your central server. You can be confident that you are satisfying transparency requirements while enhancing safety and reliability in regard to hygiene, logistics and service. The perfect complement to external monitoring activities such as routine inspections.

- A_0 value achieved
- A_0 value set
- Date and time cycle started
- Cycle duration
- Machine identification no.
- Cycle counter
- Program no.
- Temperature profile of disinfection phase
- Online remote monitoring
- Machine parameters
- and much more