

The success story continues: **TOPIC** optional extras



Customers just can't get enough of MEIKO TOPIC machines – once they have experienced these feature-packed and user-friendly systems, they find themselves wanting even more! That's why we offer such a wide array of additional features and optional extras.




- **Non-contact opening and closing – (optional)**
Automatic operation of the wash chamber door combines ease of use with maximum safety for both patients and staff.
- **Air cooling and drying (optional)**
This TOPIC upgrade uses a patented process to cool and dry care utensils with ordinary air. As well as creating a safer and more hygienic workplace, this feature can also save you up to 8 litres of water in each wash cycle. It also ensures that no steam is released into the room when you open the door at the end of the programme.
- **Disinfection of all water-carrying lines in the system (optional)**
This system uses steam to automatically disinfect the machine after each wash cycle – without incurring any extra running costs.

Please get in touch with the MEIKO team to find out how we can meet your specific requirements.

Versatile **TOPIC** configurations: perfect for any environment

TOPIC cleaning and disinfection appliances are the perfect solution for any workplace. Whether as a standalone appliance or as part of a tailor-made combined care unit – TOPIC is just the job!



MEIKO accessories

Are you making changes to your workplace environment?

Our range of genuine MEIKO accessories will help you get everything sorted! From slop sinks and rubbish bag holders to storage shelves and central drain piping, all MEIKO accessories are fully compatible with all the systems we offer.

See our MEIKO accessories brochure for more details.

Welcome to MEIKO.

