

DOYEN®

Rinse aid with water softener function and detergent for MEIKO washer-disinfectors and multiwashers



MEIKO Care



Coordinated and optimally effective

When carers reuse care utensils, one thing must be absolutely certain: The items are not only completely clean, but also perfectly clean in terms of hygiene. That is the job of our technology, along with perfectly paired chemicals.

However, there are a few things to consider when choosing the right products. What is the water like? Which stubborn residues need to be removed?

Whether DOYEN detergent or rinse aid (with water softener function) - the combination of powerful active agents provides hygiene safety in our washer-disinfectors and multiwashers, but is also gentle on materials and the environment. And that's not all: higher material protection and reduced risk of scaling reduce follow-up costs.

So what is in the canister? Most importantly: decades of experience. Experience that you can rely on.

Benefits of DOYEN products:

- Compatibility with a wide range of materials
- Specially developed for MEIKO washer-disinfectors and multiwashers
- Optimum cleaning and rinsing results
- Reduced formation of foam
- Excellent calcium binding capacity to get care utensils perfectly clean while keeping machines in great condition
- Protects the environment by minimising the amount of chemicals required



DOYEN® R100

Mildly alkaline detergent

- Powerful and intensive cleaning
- Particularly gentle on materials and eco-friendly
- Suitable for all water hardness levels

DOYEN® RS120

Acidic detergent

- Powerful and intensive cleaning
- Perfect results with ointment residues
- Suitable for all water hardness levels

DOYEN® RF90

Highly alkaline detergent

- Powerful and intensive cleaning
- For use in the multiwasher TopClean 60 (not for washer-disinfectors)
- Suitable for all water hardness levels

DOYEN® SK10

Rinse aid and water softener

- Sparkling cleaning results for care utensils
- For use in water hardness up to 10°dH
- High cost efficiency

DOYEN® SK22E

Rinse aid and water softener

- Sparkling cleaning results for care utensils
- Suitable for all MEIKO washer-disinfectors
- Can be used with all water hardnesses

DOYEN® SK22E mit Apfelduft

Rinse aid and water softener

- Sparkling cleaning results for care utensils
- Neutralises odours in the utility room
- Can be used with all water hardnesses

DOYEN® SK33E



Rinse aid and water softener

- Sparkling cleaning results for care utensils
- Suitable for demanding water conditions
- Applicable for high degrees of hardness from 30°dH or phosphate inoculations





Blue to blue, grey to grey - the colour coding and the matching moulded inserts of the suction lance and canister ensure clear use of detergent and rinse aid. There is no way to mix them up. This ensures more user and process reliability. The suction lance coding system only works with DOYEN products.



MEIKO Maschinenbau GmbH & Co. KG · Englerstr. 3 · 77652 Offenburg, Deutschland
Tel. +49 781 203 1717 · chemicals@meiko-global.com · www.meiko.com

Manufacturer:



etol Eberhard Tripp GmbH
Allerheiligenstraße 12 · 77728 Oppenau, Deutschland
Tel. +49 7804 41-127 · Fax +49 7804 41-168
www.etol.de · hygiene@etol.de

