

Discover the new pearl
in our care collection



BLUE OPERATING CONCEPT

The glass operating panel uses a progress bar, status display and our signal concept to provide an overview of machine operation at all times. Only the elements and functions lit up in blue at any moment can be used or selected.

POWERFUL CLEANING SYSTEM

The wash chamber is well lit and features a telescopic rotary jet, cleaning nozzles around 360° and a wash arm for efficient care utensil cleaning.

ADJUSTABLE A₀ VALUE OF UP TO 12,000

These washer-disinfectors can be activated with A₀ values of up to 12,000. Select the value according to your needs, such as potential outbreaks or use on isolation wards.

MEIKO SIGNAL CONCEPT

Blue = ready for operation. Green = cleaning and disinfecting. Red = important message. The door handle and lighting elements on the operating panel show care staff what the appliance's status is.

MEIKO CONNECTIVITY

The *MEIKO Connect* app allows you to view and manage all relevant data from the *MEIKO TopLine*. Stay informed even outside of the utility room through connectivity via Bluetooth, WLAN or LAN.

AUTOMATIC DOOR OPENING AND CLOSING WITH HAND SENSOR

This option allows for 100% contact-free opening and closing of the wash chamber door on the washer-disinfector. So how does it work? Simply hold your hand in front of the sensor.

LOAD DETECTION

This option suggests the correct cleaning programme for the care utensils loaded and starts it automatically after a few seconds.

AUTOMATIC DOOR OPENING AND CLOSING WITH FOOT SWITCH

A switch close to the floor for foot operation opens the wash chamber door. The ideal option to complement automatic door opening and closing with a hand sensor.

The new pearl in our care collection – technology that makes more time for people

No need to select, confirm and start the right programme – the *MEIKO TopLine* can do that. Hands full when you want to open the appliance door? No problem – it's all contact-free.* Short on time between patients? The *MEIKO TopLine* shows you from a distance whether the care utensils are ready to be unloaded.

Discover the new pearl in our care collection: the first washer-disinfector with an assistant system – this appliance boasts user convenience, intelligent functionality and future-proof hygiene settings. Clever. Tough. Equipped for the realities of everyday work in the care sector.

Discover the new pearl in the *MEIKO* care collection: www.meiko.info/pearl



* Full range of advantages only available when all options are selected



The clean solution

Member of the **MEIKO** GROUP

Manufacturer:

MEIKO Maschinenbau GmbH & Co. KG · Englerstr. 3 · 77652 Offenburg, Germany

Phone +49 781 203-0 · Fax +49 781 203-1179

www.meiko-global.com · info@meiko-global.com

