


MEIKO – TopLine Cleaning and Disinfection Appliances

MEIKO SAN 20 BW / TopLine 40



Combined sluce unit with slop sink, hand wash basin, wash basin and cleaning and disinfection appliance. Integrated with a large work surface and additional cabinet space.



Technical data

Cold water : R ½ ", 5 -25 °C
 Flow : > 18 ltr/min.
 Flow pressure : > 1 bar

Warmwasser : R ½ ", 45 – 60 °C
 Durchfluss : > 18 ltr/min.
 Fließdruck : > 1 bar

Ablauf : DN 100

Electrical installation : 3N PE ~ 400V, 50 Hz, 4.6 kW
 Alternatively:
 Electrical installation : 1N PE ~ 230V, 50 Hz, 2.8 kW

Workplace-related
 emission
 value LpA : 50 dB

We reserve the right to make technical changes as part of product improvements.