


## MEIKO – TopLine Cleaning and Disinfection Appliances

### MEIKO TopLine 40 AT (chamber door mechanism)

Cleaning and disinfection appliance as a cabinet model.  
Cabinet model fitting under a work surface.  
**Hygienic opening and closing of the door by means of a proximity switch** (pedal as an optional extra).



	<p><b>Technical data</b></p> <p>Cold water : R ½ ", 5 -25 °C Flow : &gt; 18 ltr/min. Flow pressure : &gt; 1 bar</p> <p>Hot water : R ½ ", 45 -60 °C Flow : &gt; 18 ltr/min. Flow pressure : &gt;1 bar</p> <p>Outlet pipe : DN 100</p> <p>Electrical installation : 3N PE ~ 400V, 50 Hz, 4.6 kW Alternatively: Electrical installation : 1N PE ~ 230V, 50 Hz, 2.8 kW</p> <p>Workplace-related emission value LpA : 50 dB</p>
---	---

We reserve the right to make technical changes as part of product improvements.