

TopLine product range

Cleaning and disinfection appliance





What do people all over the world expect from today's nursing profession? Certainly compliance with relevant laws and standards and a duty of care, but also proper adherence to occupational health and safety guidelines plus regular monitoring and documentation of infection control measures. The aim is clear – to prevent infectious diseases in hospitals and care homes.

MEIKO can actively assist you in this task by providing professional expertise and reliable cleaning and disinfection technology, so you get the reassurance of reliable hygiene from a single source. We can help you plan your facilities and choose the right technology and chemicals. What's more, we offer superb technical support and professional consulting services even after the sale, because we believe in forging long-term partnerships with our customers.

Whatever your needs, our bespoke hygiene management solutions have got you covered. Meeting statutory and infection control requirements is a serious business – but with MEIKO at your side you can be confident of getting it right.

Here's how MEIKO TopLine helps prevent infections and protect patients and staff:

- Disinfects the entire system every time you start a wash program
- Cools and dries care utensils with HEPA-filtered air
- Completely steam-tight system. No steam is released, not even when the door is opened at the end of the wash program
- Certified to DIN EN ISO 15883 by an independent test institute which, in turn, is accredited by the central authority of the German federal Länder for health protection regarding medicinal products and medical devices
- Energy-optimised disinfection process. The machine only uses the energy required to reach the preset A_0 value
- Displays target and actual A_0 value
- A_0 value of 60 and more can be set at will
- Helpful and informative 4 x 20-character plain text display with energy-saving mode
- Nozzle system configured to get care utensils perfectly clean, including a unique controllable telescopic rotary jet designed to thoroughly clean the inside of the utensils

MEIKO TopLine single appliances – the rapid solution for renovation and modernisation projects

Many of the dirty utility rooms of hospitals and care homes no longer comply with current regulations and hygiene standards as well as the demands of nursing staff. TopLine single appliances were specially designed for rapid and uncomplicated installation to replace

TopLine 10 – for the clear-of-the-floor wall mounting (with or without slop sink)



The easy way of installing a cleaning and disinfection appliance. Space-saving and ready for connection.

Dimensions (mm):
w 500, d 450 or 588, h (appliance) 1630



A practical and compact work unit consisting of the *TopLine 10* cleaning and disinfection appliance and a slop sink specially designed to accompany it.

Dimensions (mm):
w 1000, d 450 or 588, h (appliance) 1630

out-dated cleaning and disinfection technology. They can be modified to fit into the position of the old appliance whatever the make of the original – and to use the same connections.

TopLine 20 – floor mounted appliances (with or without slop sink)



Cleaning and disinfection appliance with plinth. Just place in position and connect. The ideal way of utilizing existing connections.

Dimensions (mm):
w 500, d 450 or 588, h (appliance) 1730



A work unit which is versatile in its functions, easy to install and compact. Consisting of the *TopLine 20* cleaning and disinfection appliance with plinth and a slop sink specially designed to accompany it.

Dimensions (mm):
w 1000, d 450 or 588, h (appliance) 1730

MEIKO TopLine single appliances – seamless and adaptable hygiene security for every work room

TopLine 30 – for the flush wall installation into bathrooms



Cleaning and disinfection appliance *TopLine 30* for the bathroom of the patients' rooms, e. g. on the isolation ward. With mounting frame.

Dimensions (mm):
Basic frame:
w 650, d 450, h 2210
Visible fronts:
w 530, h 1760



Cleaning and disinfection appliance *TopLine 30* with WC outlet for the bathroom in the patients' room, e. g. on the isolation ward. With mounting frame and complete piping as well as central connections for all commercially available WC types.

Dimensions (mm):
Basic frame:
B 650, T 450, H 2.210
Visible fronts:
w 530, h 1760

TopLine 40 – the practical table top model



The *TopLine 40* cleaning and disinfection appliance is the perfect complement to the existing cabinet system. The appliance has a large, uninterrupted work and storage surface.

Upstand available for an additional charge.

Dimensions (mm):
Appliance:
w 900, d 600, h 900



The *TopLine 40* cleaning and disinfection appliance with built-in hand wash basin and small storage surface is a suitable and small versatile combined work unit.

Upstand available for an additional charge.

Dimensions (mm):
Appliance:
w 900, d 600, h 1000
Hand wash basin:
w 340, d 370, h 150

TopLine 40 E – the fastest way to upgrade

Are you renovating or upgrading your facility? The *TopLine 40 E* is the perfect choice for quickly replacing existing appliances installed under work surfaces. Configurable to fit in all available recesses.

Dimensions (mm):
Appliance:
w 700, d 588, h 860



MEIKO TopLine combined care unit – hygiene security for all requirements

Hygiene – made to measure. We will work with you to develop an ergonomic utility room design to meet the most demanding requirements of cost-effectiveness, optimisation of work processes, health and safety and hygiene – not to mention any special factors

SAN 14 – TopLine 40 with slop sink / with slop sink and hand wash basin



Combined care unit with work and storage surface, splash screen, slop sink and upstand.

Dimensions (mm):
w 1400, d 600, h 900
Slop sink:
h 700
Sink:
Ø 362



The SAN 14 B compact combined care unit for care needs in the dirty utility room, complete with work and storage surface, splash screen, slop sink, hand wash basin and upstand.

Dimensions (mm):
w 1400, d 600, h 1000
Slop sink:
h 700
Sink:
Ø 362
Hand wash basin:
w 340, d 370, h 150

applicable to your premises. We can offer you a hygiene solution customised to meet your required hygiene standards in every room.

SAN 19 B – TopLine 40 with slop sink hand wash basin and cabinet

Combined care unit with a large work surface and a cabinet space of 500 mm. With work and storage surface, splash screen, slop sink and upstand.

Dimensions (mm):
w 1900, d 600, h 900
Slop sink:
h 700
Sink:
Ø 362
Hand wash basin:
w 340, d 370, h 150
Cabinet space:
w 500



SAN 20 BW – TopLine 40 with slop sink, hand wash basin, wash basin and cabinet

Combined care unit with small dimensions and storage surface. With 600 mm cabinet space. With work and storage surface, splash screen, slop sink and upstand.

Dimensions (mm):
w 2000, d 600, h 1000
Slop sink:
h 700
Sink:
Ø 362
Hand wash basin:
w 340, d 370, h 150
Wash basin:
w 500, d 400, h 300
Cabinet space:
w 600



MEIKO TopLine combined care units – hygiene security for all requirements

SAN 24 W – TopLine 40 with slop sink, wash basin, work surface and cabinet 1000 mm

A combined care unit with a large 1000 mm wide cabinet – a very popular choice.

Dimensions (mm):
w 2400, d 600, h 900
Slop sink:
h 700
Sink:
Ø 362
Wash basin:
w 500, d 400, h 300
Cabinet space:
w 1000



SAN 29 BW – TopLine 40 with slop sink, hand wash basin, wash basin and cabinet 1500 mm

Large combined care unit with an ideal length and basin dimensions. All cleaning tasks are optimally combined into a single unit. The basins are ergonomically arranged. The work surface is generously dimensioned.

Dimensions (mm):
w 2900, d 600, h 900
Slop sink:
h 700
Sink:
Ø 362
Hand wash basin:
w 340, d 370, h 150
Wash basin:
w 500, h 400, d 300
Cabinet space:
w 1500



MEIKO TopLine ease of operation – hands-free protection and hygiene

Ease of use and maximum hygiene security for the nursing staff as well as the greatest possible protection for your patients - all guaranteed by the hands-free versions of *TopLine* equipped with an infrared sensor to operate the door. Also available with an automatic program start on request. The efficient safety stop also protects staff against being accidentally trapped by the door. A foot switch to operate the door is a further option.

TopLine 20 & TopLine 40 – hands-free opening and closing or using the foot switch



TopLine 20 floor model and *TopLine 40* table top model with automatic door-opening (optional) and foot switch (optional).



Opening and closing by means of a proximity switch (optional).



Opening and closing by using the foot switch (optional).

Accessories and range of fixtures

In addition to the MEIKO cleaning and disinfection technology, we can also supply you with a comprehensive range of accessories for safe and ergonomic working and the hygienic storage of care utensils. We illustrate below 6 examples of practical solutions. We would be happy to develop a range of suggested solutions specially for you. Just give us a call. We would be pleased to help you at a moment's notice.



Drying and storage racks*
for disinfected urine bottles, bedpans and receptacles, with storage grid or drip tray. Available in widths of 500 mm to 1400 mm.



Cupboards with doors**
Cupboard units with hinged doors and one intermediate shelf. Available in widths of 500 mm, 600 mm, 800 mm, 1000 mm, 1200 mm and 1500 mm.



Wall shelving*
Rack with height-adjustable shelves available in all shelf lengths. Shelf depth: 350 mm.



Cupboards without doors
Cupboards for easy storage with one intermediate shelf. Available in the following widths: 500 mm, 1000 mm and 1200 mm.



Cleaning combination*
Cupboards for installation under work surfaces (complete with drawers and other facilities) can be designed in a range of widths and heights for combined care units.



Slop sink*
Standard model available as an option with integral rim flushing which can be operated independently of the cleaning and disinfection appliance.

* Pictures contain special equipment
** Optional: lamp shade and lamp

